



Bottle Label - BONAVERA COUNT EZ
Ver 1.0 Date 12.06.2018
SIZE 120MM x 65MM
MFG BY VECTOR BIOTEK PVT. LTD.

 <p>BONAVERA Count EZ</p>	<p>BONAVERA COUNT EZ</p>	<p>Code : BCC01 100ML</p>
<p>Precautions:</p> <ol style="list-style-type: none">1. Do not inhale the reagent or pipette by mouth.2. Avoid contact with skin/eyes. In case of contact, wash off immediately with plenty of water and seek medical attention as early as possible.3. Do not use if package is damaged.4. Dispose waste product in accordance with local Govt. regulations.5. Cap the bottle tightly after use.	<p>Lot No. : Mfg. : Exp. : MRP Rs. : (Inclusive of All Taxes)</p>	<p>Storage and Stability: The unopened reagent is stable up to the expiry date mentioned on the label when stored between 15°C - 30°C. Use BONAVERA COUNT EZ within 6 months of first opening.</p>
<p>Store at R.T. For Laboratory Use Only Not for Medicinal Purpose</p>	<p>Marketed By:</p>	<p>Directions of Use: Open the container cover of BONAVERA COUNT EZ during the system shutdown, press the aspirate key to aspirate BONAVERA COUNT EZ, as per instructions given in the user manual of the Analyzer.</p>
	<p>Plot No.12 Survey No. 169 (Sibdrn) Near Plot No. 425, New GIDC, Kabilpore, Navsari - 396 424, (GUJ), India Ph : +91-93 757 94111 Fax : +91-26 37-26 5105 E-mail : biogeny@beaconindia.com Web : www.biogenydiagnostics.com</p>	<p>Intended Use: BONAVERA COUNT EZ is used for cleaning the tubes & bath before shutdown in BONAVERA COUNT 3-part differential hematology analyzers.</p>
		<p>Active Ingredients: Sodium Chloride <0.85% Sodium Citrate di-hydrate <0.25% Calcium Chloride di-hydrate <0.03% Enzymes <0.65% Dye <0.01%</p>
		<p><small>VECTOR BIOTEK PVT. LTD. 425, NEW GIDC, KABILPORE, NAVSARI - 396424 (GUJ) INDIA</small></p>